



# EC Type-Examination Certificate

**No. CE 572204**

*The British Standards Institution hereby grants to:*

**3M United Kingdom Plc  
3M Centre  
Bracknell  
Berkshire  
RG12 8HT  
United Kingdom**

*In respect of:*

**Respiratory protective devices to EN 149:2001+A1:2009  
Filtering half masks to protect against particles  
Model: 3M™ Aura™ 1883+ Particulate Respirator.**

On the basis of our examination, under the requirements of Council Directive 89/686/EEC "Personal Protective Equipment" Article 10, EC Type-Examination.

*For and on behalf of BSI, a Notified Body for the above Directive (Notified Body Number 0086):*

*Rob Wallis, Chair Certification Body Management Committee*

First Issued: **28/10/2011**

Latest Issue: **28/10/2011**

*raising standards worldwide*™

Page 1 of 3



This certificate remains the property of The British Standards Institution and shall be returned immediately upon request. To check its validity telephone +44 (0)845 080 9000

The British Standards Institution, Kitemark Court, Davy Avenue, Knowlhill, Milton Keynes MK5 8PP United Kingdom. Tel: +44 (0)845 080 9000, Web: [www.bsigroup.com/certification](http://www.bsigroup.com/certification). The British Standards Institution Headquarters: 389 Chiswick High Road, London W4 4AL Tel: +44 (0)20 8998 9001. The British Standards Institution is incorporated by Royal Charter.



## No. CE 572204

### Product Details

<b>Product name:</b>	<b>Particulate Respirator.</b>
<b>Model:</b>	<b>3M™ Aura™1883+ (Health Care market).</b>
<b>Product type:</b>	<b>Filtering half mask to protect against particles.</b>
<b>Technical Specification:</b>	<b>Harmonized European Standard EN 149:2001+A1:2009.</b>
<b>EN 149 classification:</b>	<b>FFP3. Single shift use – symbol NR. Dolomite clogging option – symbol D.</b>
<b>Additional protection characteristics:</b>	
<b>Technical Specification:</b>	<b>Harmonized European Standard EN 14683:2005.</b>
<b>EN 14683 classification:</b>	<b>Type IIR.</b>

The product assessments were based on the following Standards:  
BS EN 149:2001+A1:2009, the English language version of EN 149:2001+A1:2009, respiratory protective devices – filtering half masks to protect against particles, both documents incorporating Corrigendum dated July 2002;  
and BS EN 14683:2005, the English language version of EN 14683:2005.

The shrouded valved respirator is designed to protect against solid and non-volatile liquid particles, is held onto the face by a pair of elasticated straps and is a single shift device.

3M has Self Certified this product to the requirements of Annex VII of the European Community Directive 93/42/EEC (Medical Devices Directive) as a Class 1 device.



## No. CE 572204

### Certification Administration Details

Technical File Reference: TF0266


### Certificate Amendment Record and internal BSI report related to this Certificate:

Issue date	Comments	BSI Report No.
See 'Latest Issue' below	First issue.	0086:7650930

### Certificate validity

The Certificate holder is responsible for ensuring that the Notified Body is advised of changes to any aspect of the overall processes utilised in the manufacture of the product, failure to do so could invalidate the Certificate in respect of product manufactured following the introduction of such changes.

The validity of the Certificate is also dependent on the maintenance of the EC quality of production by monitoring system, PPE Article 11B, as referenced on BSI Certificate 67902.



First Issued: **28/10/2011**

Latest Issue: **28/10/2011**

Page 3 of 3

*raising standards worldwide*<sup>TM</sup>



This certificate remains the property of The British Standards Institution and shall be returned immediately upon request. To check its validity telephone +44 (0)845 080 9000

The British Standards Institution, Kitemark Court, Davy Avenue, Knowlhill, Milton Keynes MK5 8PP United Kingdom. Tel: +44 (0)845 080 9000, Web: [www.bsigroup.com/certification](http://www.bsigroup.com/certification). The British Standards Institution Headquarters: 389 Chiswick High Road, London W4 4AL Tel: +44 (0)20 8998 9001. The British Standards Institution is incorporated by Royal Charter.